

## MANUAL:

### **Cas 13 Activity Assay Protocol**

Catalog Number: orb1742853

#### **Reaction Components:**

<b>Active Cas13</b>	It is recommended that the researcher perform a serial dilution of Active Cas13 diluted with 1X Reaction Buffer for optimal results.
<b>Cas13 Reaction Buffer (10x) (Cat #: orb1742843)</b>	Buffer components: 200 mM HEPES pH 7.0, 500 mM KCl, 50 mM MgCl <sub>2</sub> , 1 mg/mL BSA.
<b>CRISPR RNA</b>	Single guide RNA containing a target-specific spacer sequence and a scaffold sequence responsible for recruiting the Cas enzyme.
<b>Target RNA activator</b>	An RNA oligo that contains the complementary sequence to the crRNA.
<b>Reporter Substrate:</b>	A synthetic RNA oligonucleotide labelled with a fluorophore and a quencher on the two termini.

#### **Assay Protocol:**

The RNA trans-cleavage activity of Cas13 is detected in a CRISPR-based fluorescent reporter assay. RNA-guided RNA binding to Cas13 activates the enzyme, inducing target RNA cutting as well as collateral RNase activity. The latter leads to the degradation of the reporter substrate that emits a light signal on cleavage.

##### **Step 1.**

Thaw the active Cas13 on ice. Prepare 1X Reaction Buffer. Equilibrate the buffer, guide RNA, target RNA, and reporter substrate to ambient temperature.

##### **Step 2:**

Prepare the following working solutions with 1X Reaction Buffer:

- 4X final concentration of Active Cas13

- 4X final concentration of guide RNA
- 2X final concentration of activator/ reporter mix containing the RNA activator and reporter substrate.

### **Step 3:**

In a half-area solid black 96-well plate, add the following components and pre-incubate at room temperature for 10 minutes:

- Component 1. 10  $\mu$ L of 4X Active Cas13
- Component 2. 10  $\mu$ L of 4X guide RNA

*Note: A blank control can be set up as outlined in step 3 by replacing the enzyme working solution with an equal volume of the assay buffer*

### **Step 4:**

To each assay well, add 20  $\mu$ L of the 2X activator/reporter mix. Shake the plate for 1 minute on a tabletop orbital shaker. Seal the assay wells with microplate sealing tape and incubate at 37°C for 15-30 minutes.

### **Step 5:**

Equilibrate the plate to ambient temperature and then remove the microplate sealing tape. Read fluorescence on a microplate reader.