

FlysisAmp Cells Lysis Kit

CL101



Instruction for Use

Version 24.1

Contents

01/Product Description	02
02/Components	02
03/Storage	02
04/Applications	02
05/Self-prepared Materials	02
06/Notes	02
07/Mechanism & Workflow	04
08/Experiment Process	05
09/FAQ & Troubleshooting	07

01/Product Description

The FysisAmp Cells Lysis Kit is a fast, efficient kit for harvesting cellular RNA and is compatible with $10^5 - 10^6$ cultured cells (adherent or suspension cells). Genomic DNA elimination and RNA acquisition can be completed directly in the cell culture plate, and the whole process only takes 7 min. This method does not require the steps of heating, washing, tube transfer, and centrifugation in traditional extraction methods, which improves efficiency while also reducing losses during sample purification. This kit is compatible with a wide range of cell types and cell quantities, and can also solve the problem of unsatisfactory extraction results with low cell numbers in traditional extraction methods. This kit can be used for high-throughput (96-/384-well cell culture plate) acquisition of cellular RNA and for various downstream experiments such as RT-PCR and RT-qPCR.

02/Components

Components		CL101-01 100 rxns (50 μ l/rxn)	CL101-02 500 rxns (50 μ l/rxn)
<input type="checkbox"/>	FysisAmp Cells Lysis Buffer*	5 ml	25 ml
<input checked="" type="checkbox"/>	DNase I	200 μ l	1 ml
<input checked="" type="checkbox"/>	Enhancer Solution	200 μ l	1 ml
<input checked="" type="checkbox"/>	FysisAmp Cells Stop Buffer	500 μ l	2 \times 1.25 ml

* FysisAmp Cells Lysis Buffer can be stored stably for 1 year at $2 \sim 8^\circ\text{C}$ after thawing.

03/Storage

Store at $-30 \sim -15^\circ\text{C}$ and transport at $\leq 0^\circ\text{C}$.

04/Applications

It is applicable for RNA extraction from $10^5 - 10^6$ cultured cells (adherent or suspension cells).

05/Self-prepared Materials

1 \times PBS buffer, 1.5 ml RNase-free centrifuge tube, RNase-free pipette tip, 0.2 ml RNase-free eight-tube PCR strip/PCR tube, etc.

06/Notes

For research use only. Not for use in diagnostic procedures.

1. Before use, please check whether there is any ice in the reagent. If any component freezes, place it at room temperature to dissolve, and mix well before use.

2. If there are too many cells during the lysis process, it may lead to incomplete lysis or inhibit subsequent reverse transcription and qPCR amplification. Therefore, before operation, adherent cells need to select the corresponding volume of lysis working solution according to the number of plate wells; suspension cells need to be estimated number of cells to select the corresponding volume of lysis working solution.
3. When using fresh samples, if they cannot be lysed in time, wash the cells once with pre-cooled PBS and store them at $-85 \sim -65^{\circ}\text{C}$, and avoid repeated freezing and thawing. To avoid RNA degradation, cells should be processed and stored as quickly as possible.
4. Gently pipette 8 - 10 times to mix well. Avoid air bubbles caused by vigorous shaking.
5. When using this kit, please wear lab coats, disposable latex gloves, disposable masks, and use RNase-free consumables to avoid RNase contamination.
6. The lysis process is performed in an RNase-free environment; unless otherwise specified, operations are performed at room temperature ($15 \sim 25^{\circ}\text{C}$).
7. A microplate shaker can be used to shake at 600 rpm for 30 sec instead of pipetting to mix.

07/Mechanism & Workflow




Fig 1. Workflow of FysisAmp Cells Lysis Kit

08/Experiment Process

◇ Adherent cells:

1. Sample processing and preparation

Gently remove the cell culture medium from the adherent cells in different well plates, then slowly add an equal volume of pre-cooled 1 × PBS to the plate to wash the adherent cells. Remove the PBS carefully, and place the cells on ice for later use.

▲ When removing culture medium or PBS, the action should be as gentle as possible to avoid cell loss during the washing process.

2. RNA acquisition

a. Prepare lysis working solution (take 96-well plate as an example)

Components	Volume	
FlysisAmp Cells Lysis Buffer	46 μl	<input type="checkbox"/>
DNase I	2 μl	<input checked="" type="checkbox"/>
Enhancer Solution	2 μl	<input checked="" type="checkbox"/>
Total	50 μl	

▲ Lysis working solution should be prepared according to the required volume of the well plate, the ratio of each component is FlysisAmp Cells Lysis Buffer : DNase I : Enhancer Solution = 23 : 1 : 1.

▲ After preparation, mix by inverting 10 - 15 times, avoid violent vortexing; store on ice after mixing, please use within 1 h.

b. Refer to the table below and add the corresponding lysis working solution to the aliquoted cell sample, pipette gently 8 - 10 times to mix thoroughly, and let stand at room temperature for 5 min to lyse the cells.

Components	Volume					
Cell culture plate	384 well	96 well	48 well	24 well	12 well	6 well
Lysis working solution volume	25 μl	50 μl	100 μl	150 μl	250 μl	750 μl

c. Add FlysisAmp Cells Stop Buffer as shown in the table, pipette gently 8 - 10 times to mix, and let stand at room temperature for 2 min to terminate the reaction.

Components	Volume					
Cell culture plate	384 well	96 well	48 well	24 well	12 well	6 well
FlysisAmp Cells Stop Buffer	2.5 μl	5 μl	10 μl	15 μl	25 μl	75 μl
Total	27.5 μl	55 μl	110 μl	165 μl	275 μl	825 μl

▲ The lysate can be used for subsequent detection experiments, such as reverse transcription, one-step RT-qPCR/RT-PCR.

▲ The lysate can be stored on ice for 2 h. If long-term storage is required, it should be stored at -80°C.

◇ Suspension cells:

1. Sample processing and preparation

a. Remove the culture medium: Transfer the suspension cells to a centrifuge tube, centrifuge at 4°C, 1,000 rpm (930 × g) for 5 min, collect to the bottom of the tube, and remove the culture medium.

- b. Wash cells: Add pre-cooled 1 × PBS (500 µl/10⁵ cells), pipette gently 8 - 10 times until there are no obvious cell clusters; centrifuge at 1,000 rpm (930 × g) at 4°C for 5 min. Collect to the bottom of the tube and remove PBS.
- c. Dilute cells: According to the experimental requirements, dilute the cells with pre-cooled 1 × PBS, distribute the diluted cell suspension into eight-tube strips, 5 µl per tube (keep the number of cells at 10 - 10⁵ cells/tube), and place on ice.
- ▲ When removing culture medium or PBS, the action should be as gentle as possible to avoid cell loss during the washing process.

2. RNA acquisition

a. Preparation of lysis working solution

Components	Volume		
	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
Cell number	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
FlysisAmp Cells Lysis Buffer	24 µl	46 µl	184 µl
DNase I	1 µl	2 µl	8 µl
Enhancer Solution	-	2 µl	8 µl
Total	25 µl	50 µl	200 µl

▲ After preparation, mix by inverting 10 - 15 times, avoid violent vortexing; store on ice after mixing, please use within 1 h.

- b. Refer to the table below and add the corresponding lysis working solution to the aliquoted cell sample, pipette gently 8 - 10 times to mix thoroughly, and let stand at room temperature for 5 min to lyse the cells.

Components	Volume		
	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
Cell number	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
Lysis working solution	25 µl	50 µl	200 µl
Total	30 µl	55 µl	205 µl

- c. Add FlysisAmp Cells Stop Buffer as shown in the table, pipette gently 8 - 10 times to mix, and let stand at room temperature for 2 min to terminate the reaction.

Components	Volume		
	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
Cell number	10 - 10 ³	10 ³ - 10 ⁵	10 ⁵ - 10 ⁶
FlysisAmp Cells Stop Buffer	2.5 µl	5 µl	20 µl
Total	32.5 µl	60 µl	225 µl

▲ The lysate can be used for subsequent detection experiments, such as reverse transcription, one-step RT-qPCR/RT-PCR.

▲ The lysate can be stored on ice for 2 h. If long-term storage is required, it should be stored at -80°C.

09/FAQ & Troubleshooting

FAQ	Reasons	Solutions
No PCR product or low yield	1. Insufficient cell lysis	1. Insufficient addition of Lysis Buffer: Follow the recommended number of cells and well plates in the instruction. 2. Uneven mixing of cells and Lysis Buffer: When adding Lysis Buffer, extend the pipette tip below the liquid level. Mix by pipetting up and down 8 - 10 times.
	2. Stop Buffer is not added or mixed evenly	When adding Stop Buffer, extend the pipette tip below the liquid level. Mix by pipetting up and down 8 - 10 times.
	3. Insufficient cell input or cell lysate input	Increase the amount of cell input or increase the amount of lysate input according to the instructions.
	4. Excessive cell input: Excessive cell contents inhibit amplification.	Perform the experiment according to the instruction.
	5. The sample does not contain the target RNA or the expression level of target RNA is low	Increase the amount of cell input or adjust the amplification system.
	6. The frozen cells were rinsed with PBS	After freezing, the cells have broken. PBS washing results in RNA loss. It is recommended to re-culture the cells and repeat the experiment.
RNA degradation	1. RNA was degraded before the experiment	Place fresh cells on ice promptly after processing to avoid prolonged placement at room temperature.
	2. The lysate was left at room temperature for too long	The lysate should not be allowed to stand at room temperature for more than 20 min or on ice for more than 2 h. Follow-up tests should be carried out in a timely manner or the lysate should be frozen at -80°C.
	3. Excessive cell input: Cellular RNases degrade RNA	Reduce the amount of cell input and perform the experiment according to the instruction.
Genomic DNA residues	1. Excessive cell input	Cell input should be $<10^6$.
	2. Insufficient cell lysis	Extend the lysis time appropriately, but not more than 10 min.



Nanjing Vazyme Biotech Co.,Ltd.

Tel: +86 25-83772625

Email: info.biotech@vazyme.com

Web: www.vazyme.com

Loc: Red Maple Hi-tech Industry Park, Nanjing, PRC

Follow Us

