

Product Datasheet

2 × Rapid Taq Master Mix (P222-02)

Description	2 × Rapid Taq Master Mix contains Taq DNA Polymerase, elongation promoting factor, dNTP and an opti...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. It can be stored at 4°C for 3 months after thawing.
Note	For research use only
Application notes	Rapid amplification: 15 sec/kb; extension speed up to 1 sec/kb (≤ 1kb fragments) Convenience: Include tracking dye, direct loading of PCR product on gel Good stability: Stable performance after 50 freeze-thaw cycles Ready-to-use: Simplified workflow and minimized pipetting steps with 2 × master mix format
Expiration Date	12 months from date of receipt.

