


Product Datasheet

SuperKine West Pico PLUS Chemiluminescent Substrate (orb1147873)

Description	ECL
Storage	Stable for one year at 4-8°C from date of shipment. Exposure to the sun or any other intense light can harm the working solution. Short-term exposure to typical laboratory lighting will not harm the working solution.
Note	For research use only
Application notes	SuperKine West Pico PLUS Chemiluminescent Substrate provides a highly sensitive color developing solution to enhance the chemical signal of immunological experiments. In order to obtain the best experimental results, you need to optimize all your experimental factors, including the number of samples, antibody concentration, the use of membranes and blocking agents.
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.



SuperKine™
West Pico PLUS
Chemilumines...