



Product Datasheet

Lymphoid Tissue Tissue Slide (Lymph Node -Normal) [Out of stock] (orb1227301)

Description Lymphoid Tissue Tissue Slide

(Lymph Node - Normal) [Out of

stock]

Tested IHC

Applications

Storage

Stable for up to two years when

stored in a cool, dry lab area.

Note For research use only

Application notes Blorbyt's human tissue slides

> can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.In order to avoid false negative results from DNA or RNA sequences not being present in a 4 µm cellular cross-

section, 10 µm tissue sections are recommended. 10 µm single tissue slides are available upon request. This product has not been evaluated for effectiveness in ISH or FISH applications.

Expiration Date 12 months from date of receipt.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom