

Product Datasheet

Pancreas Tissue Slide (Normal) [Out of stock] (orb1227328)

| | |
|----------------------------|--|
| Description | Pancreas Tissue Slide (Normal) [Out of stock] |
| Tested Applications | IHC |
| Storage | Stable for up to two years when stored in a cool, dry lab area. |
| Note | For research use only |
| Application notes | Biorbyt's human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol. In order to avoid false negative results from DNA or RNA sequences not being present in a 4 µm cellular cross-section, 10 µm tissue sections are recommended. 10 µm single tissue slides are available upon request. This product has not been evaluated for effectiveness in ISH or FISH applications. |
| Expiration Date | 12 months from date of receipt. |