

Product Datasheet

CXCL1 Antibody (orb1272521)



www.biorbyt.com

Describitionnts.

CXCL1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, NeA, WB

Applications

Immunogen

Produced from sera of rabbits pre-immunized with highly

pure (>98%) recombinant hGRO/MGSA (human

GRO/MGSA).

Target

CXCL1

Form/Appearance

Lyophilized

Concentration

batch dependent

Storage

GRO/MGSA antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid

repeated freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Uniprot ID

P09341

NCBI

P09341

Dilution Range

Neutralization: To yield one-half maximal inhibition [ND50] of the biological activity of hGRO/MGSA 100 ng/mL), a concentration of 1.5 - 3.0 µg/mL of this antibody is required. ELISA:To detect hGRO/MGSA by direct ELISA (using 100 µL/well antibody solution) a concentration of at least $0.5 \mu g/mL$ of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hGRO/MGSA. Sandwich: To detect hGRO/MGSA by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with our Biotinylated Anti-Human GRO/MGSAas a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hGRO/MGSA. Western Blot:To detect hGRO/MGSA by Western Blot analysis this antibody can be used at a concentration of $0.1 - 0.2 \mu g/mL$. Used in conjunction with compatible secondary reagents the detection limit for recombinant hGRO/MGSA is 1.5 - 3.0 ng/lane, under either reducing or

To detect Human GRO/MGSA by sandwich ELI...



To detect Human GRO/MGSA by Western Blot...



To detect Human GRO/MGSA by Western Blot...

Biorbyt LLC.

non-reducing conditions