

## Product Datasheet

### Antifade Mounting Medium (orb1473253)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Antifade Mounting Media   |
| <b>Tested Applications</b> | IF  |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Protect from light. Avoid freeze/thaw cycles.  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | <p>1. Bring the vial to room temperature. 2. Rinse slide to be mounted with distilled or deionized water, touch the edges of slide on a paper towel to remove excess water. Place slides on a flat surface. 3. Apply 3-4 drops of mounting medium directly on top of the specimen and spread out evenly by tilting slide back and forth or spread evenly with a 0.2 ml plastic pipette tip making sure the tissue is not touched. Excess medium can be removed by touching the edges of slide against paper towel. 4. Let stand at room temperature for about 5 minutes. 5. Apply cover slip carefully avoiding air bubbles. 6. The specimen is ready for visualization under a microscope.</p> |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |