

Product Datasheet

Human Progesterone Associated Endometrial Protein (PAEP) ELISA Kit (orb1654818)

Description Human Progesterone Associated Endometrial Protein (PAEP) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 4.69-300 ng/mL. Sensitivity: 1.9 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.