

Product Datasheet

Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit (orb1663220)

Description	Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.49 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.