

Product Datasheet

Rabbit Atrial Natriuretic Peptide (ANP) ELISA Kit (orb1667883)

Description	Rabbit Atrial Natriuretic Peptide (ANP) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.