

Product Datasheet

Uterus-Corpus Membrane Lysate (orb1672356)

Description	Uterus-Corpus Membrane Lysate
Tested Applications	WB
Preservatives	HEPES pH 7.9, MgCl ₂ , KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, and NP-40.
Form/Appearance	Liquid
Concentration	2 mg/mL or better
Storage	Store at 2-8°C for continuous use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt.
Note	For research use only
Application notes	Uterus: Corpus membrane protein lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel.
Expiration Date	12 months from date of receipt.