

Product Datasheet

Yeast DNA Preparation - Column Kit (orb1733627)

Description	The spin column based Yeast DNA Preparation Kit is designed for rapid and high purity isolation of g...
Storage	store at ambient temperature
Note	For research use only
Application notes	The obtained DNA is suitable for a variety of applications, including real-time PCR, southern blot analysis, genotyping and discovery or validation of SNP/SSR markers. Before start, add the following components as indicated on the re-spective bottle/tube: Lyticase Buffer to Lyticase. double distilled water to RNase A and Proteinase K. 96-99% Ethanol to the Washing Buffer.
Expiration Date	12 months from date of receipt.