



Product Datasheet

Mouse COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit (orb1949575)

biorbyt

Descriptionnts.

Mouse COMP(Cartilage Oligomeric Matrix

www.biorbyt.com

Reactivity	Mouse	
Range	1.57-100ng/mL	· • /
Target	СОМР	1 J J J J J J J J J J J J J J J J J J J
Note	For research use only	
Application notes	This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse COMP. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse COMP and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse COMP, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm \pm 2 nm. The OD value is proportional to the concentration of Mouse COMP. You can calculate the concentration of Mouse COMP in the samples by comparing the OD of the samples to the standard curve.	
Sample Types	serum, plasma, Tissue homogenate and Other biological samples	
Uniprot ID	Q9R0G6	
Sensitivity	0.94 ng/mL	
Expiration Date	Please enquire.	

 \mathbf{F}_{i}

Biorbyt Ltd

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom Email: **info@biorbyt.com** | Phone: **+44 (0) 1223 859-353** | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States Email: **info@biorbyt.com** | Phone: **+1 (415) 906-5211** | Fax: +1 (415) 651-8558