

Product Datasheet

A431 cells stimulated with EGF (orb198380)

Description	A431 cells stimulated with EGF
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Form/Appearance	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Suggested use: A volume of 10-25ul is suggested for immunoblotting.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.