



Product Datasheet

CO1A2 antibody (orb345832)



www.biorbyt.com

Descriptionnts. CO1A2 antibody

Species/Host Rabbit

Reactivity Bovine, Human

Conjugation FITC

Tested ELISA, FC, IF, IHC, IP, WB

Applications

Anti-Collagen Type II antibody was produced by **Immunogen**

> repeated immunizations with Collagen Type II from adult human knee cartilage and bovine

nasal cartilage.

Preservatives 0.01% (w/v) Sodium Azide

Lyophilized Form/Appearance

Concentration 1.0 mg/ml

Storage Store vial at 4° C prior to restoration. Restore

> with 0.05 mL of deionized water (or equivalent). For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room

> temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one

(1) year from date of restoration.

Note For research use only

Application notes Anti-Collagen Type II Fluorescein Conjugated

> Antibody is suitable for western blotting, IP, IHC and for ELISA. Researchers should determine optimal titers for applications that are not

stated below.

Isotype IgG

Clonality Polyclonal

Purity This product has been prepared by

> immunoaffinity chromatography using immobilized antigens followed by extensive cross-adsorption against other collagens. Some class specific anti-collagens may be specific for three-dimensional epitopes which may result in diminished reactivity with denatured collagen or formalin-fixed, paraffin embedded tissues. This antibody reacts with most mammalian Type II collagens and has expected cross-reactivity with Tyne IV and VI and negligible cross

Biorbyt Ltd.

Biorbyt LLC.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

68 TW Alexander Drive
br>Research Triangle Park
br>Durham, North Carolina < br > 27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558