

## Product Datasheet

### Facs And If Staining Buffer With Fbs [Out of stock] (orb348594)

<b>Description</b>	Facs And If Staining Buffer With Fbs
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC
<b>Preservatives</b>	0.1% (w/v) Sodium Azide
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1X
<b>Storage</b>	Store container at 4° C prior to opening.
<b>Note</b>	For research use only
<b>Application notes</b>	Fluorescent Activated Cell Sorting (FACS) and Immunofluorescence (IF) staining buffer is adequate for cell staining, washing and dilution in flow cytometry and immunofluorescence experiments.
<b>Purity</b>	FACS and IF Staining Buffer was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
<b>Hazard Information</b>	Non-Toxic
<b>Dilution Range</b>	FC: User Optimized, IHC: User Optimized, IF: User Optimized
<b>Expiration Date</b>	12 months from date of receipt.