

Product Datasheet

Facs And If Staining Buffer With Bsa (orb348598)

Description	Facs And If Staining Buffer With Bsa
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Preservatives	0.1% (w/v) Sodium Azide
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at 4° C prior to opening.
Note	For research use only
Application notes	Fluorescent Activated Cell Sorting (FACS) and Immunofluorescence (IF) staining buffer is adequate for cell staining, washing and dilution in flow cytometry and immunofluorescence experiments.
Purity	FACS and IF Staining Buffer was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Dilution Range	FC: User Optimized, IHC: User Optimized, IF: User Optimized
Expiration Date	12 months from date of receipt.