

www.biorbyt.com

Product Datasheet

A-172 Whole Cell Lysate (orb348708)

biorbyt

www.biorbyt.com

| Description ^{nts.} | A-172 Whole Cell Lysate | |
|-----------------------------|--|--|
| Conjugation | Unconjugated | |
| Tested Applications | ChIP, IP, SDS-PAGE, WB | |
| Preservatives | Preservative: None. Stabilizer: None. 1X RIPA Buffer with HALT Protease and Phosphatase Inhibitors | |
| Form/Appearance | Liquid (sterile filtered) | |
| Concentration | 4.0 mg/ml | |
| Storage | Store vial at -70° C or COLDER. For extended storage, aliquot lysate to minimize freeze/thaw cycles. | |
| Note | For research use only | |
| Application notes | orb348708 has been tested by SDS-PAGE. Multi-purpose A-172 Whole Cell Lysate are especially prepared as positive control for multiple assays including western blot, immunoprecipitation (IP), capture ELISA or other assays requiring native protein sample. For separation by SDS-PAGE and subsequent western blot analysis, lysates should be diluted by user to desired concentration in SDS-PAGE buffer with 2-mercaptoethanol or dithiothreitol as the reducing agent and heated to 95° C for 5 minutes. Sample is ready for use in immunoprecipitation and ELISA experiments, conditions should be optimized by the user. Rockland recommends its TrueBlot IP reagents for immunoprecipitation experiments. | |
| Purity | A-172 cells were grown in Dulbecco's medium supplemented with 10% fetal bovine serum. Cells were washed with PBS and then incubated on ice in modified RIPA buffer, containing 150 mM sodium chloride, 50 mM Tris HCl, pH 7.4, 1 mM EDTA, 1.0% NP-40, 0.5% sodium deoxycholic acid , 0.1% SDS and 0.01% (w/v) sodium azide to lyse the cells. Protein integrity was ensured using a cocktail of protease inhibitors with broad specificity for the inhibition of aspartic, cysteine, and serine proteases as well as aminopeptidases (0.1 mM AEBSF HCl, 0.08 µM Aprotinin, 5 µM Bestatin, 1.5 µM E-64, 2 µM Leupeptin Hemisulfate, 1 µM Pepstatin | |

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom Email: **info@biorbyt.com** | Phone: **+44 (0) 1223 859-353** | Fax: +44 (0)1223 280 240 Biorbyt LLC. 68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

mg/mi with 1X KIPA buller.