

Product Datasheet

**hsa-mir-1271-5p Real-Time RT-PCR Detection and
cel-mir-39-3p Calibration Kit [Out of stock]
(orb390601)**

Description

Detect mature miRNA hsa-mir-1271-5p

Tested Applications

PCR

Target

hsa-mir-1271-5p Real-Time RT-PCR
Detection and cel-mir-39-3p Calibration Kit

Storage

Shipped at 4°C. Store at -20°C for 6 months.

Note

For research use only

Application notes

The Real-Time RT-PCR detection and cel-mir-39-3p calibration Kit contains the whole reagent including the stem-loop RT primer & specific PCR primer set for microRNA quantification and cel-mir-39-3p primer set for calibration. It provides a sensitive and specific method of real-time PCR for relative quantification of miRNAs from total RNA samples. It is ideal for validating microRNA expression profiling data and high-throughput applications. Trait lies in even highly homogenized miRNAs can be accurately discriminated ensured by stemloop RT primer and high specific miRNAs primer set.

Experiment Notes: PCR Method: 2-Step qRT-PCR

Detection Range: Seven orders of magnitude

Sample Types: miRNA, Total RNA

Platform: ABI 7500 System/ViiA™ 7 System/StepOne™/StepOnePlus™/7900HT System; BioRad IQ5.0/CFX Connect; Stratagen MX3000p/4000p; Roche LightCycler 480

Components:
50 rxns: It contains all reagents required to perform 50 cDNA reactions and 50 PCR reactions.
100 rxns: It contains all reagents required to perform 100 cDNA reactions and 100 PCR reactions.
200 rxns: It contains all reagents required to perform 200 cDNA reactions and 200 PCR reactions.

Sensitivity

10 copies of a standard mimics and as little as 25 pg of total RNA

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558