

Product Datasheet

Bovine PIINP ELISA Kit (orb441568)

Description	Bovine Procollagen II N-Terminal Propeptide (PIINP) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	34pg/mL
Expiration Date	Please enquire.