


Product Datasheet

BlueElf Prestained Protein Marker (orb533732)

Description	The BlueElf Prestained Protein Marker is a three-color protein standard with 13 pre-stained purified...
Form/Appearance	liquid (supplied in 20 mM Tris-phosphate pH 7.5 , 2% SDS, 0.2 mM DTT, 3.6 M urea and 15% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Visualization of protein migration during SDS-polyacrylamide gel electrophoresis. Verification of protein transfer onto membranes after Western blotting. Approximate sizing of proteins on SDS-polyacrylamide gels and Western blots. Recommended Load: 3-5 µl Protein Molecular Weight Marker per lane.
Purity	> 95% (SDS-PAGE)
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



SDS-Page
analysis of
BlueElf
Prestained

...