

Product Datasheet

Equine TIMP2 ELISA Kit (Ready to Use) (orb549498)

Description	Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Equine
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	077717
Sensitivity	1.49ng/mL
Entrez	100034134
Expiration Date	Please enquire.