

Product Datasheet

Other DHEA ELISA Kit (Ready to Use) (orb549545)

Description	General Dehydroepiandrosterone (DHEA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	All
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.023ng/mL
Expiration Date	Please enquire.