



Product Datasheet

BirA antibody (orb750518)



www.biorbyt.com

Descriptionnts. BirA antibody

Species/Host Rabbit

Conjugation Unconjugated

Tested
Applications

ELISA, WB

Immunogen Anti-Biotin Ligase Epitope Tag antibody was

prepared from whole rabbit serum produced by repeated immunizations with a synthetic avi-tag peptide GLNDIFEAQKIEWH (minimal BirA

substrate; biotin is added to the sequence at the

lysine).

Preservatives 0.01% (w/v) Sodium Azide

Form/Appearance Liquid (sterile filtered)

Concentration 70 mg/ml

Storage Store vial at -20° C prior to opening. Aliquot

contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate

use.

Note For research use only

Application notes Anti-Biotin Ligase Epitope Tag Antibody has been

tested in ELISA and Western Blot. Western blot analysis of B. subtilis strains producing either (BAP-tagged) SPOIIIAH (XWB49) or (BAP-tagged) SpoIIQ (XWB46) showed strong and specific reactivity to the antisera. Minor background bands of unknown identity may be observed. The

antisera will recognize the BAP/BLT tag

irrespective of the whether the tag is biotinylated

(data not shown).

Isotype Antiserum

Clonality Polyclonal

Purity This antiserum detects the minimal peptide

substrate (GLNDIFEAQKIEWH) for the E. coli biotin ligase tag (BLT of BirA). Western blot analysis of B. subtilis strains producing either (BAP-tagged) SPOIIIAH (XWB49) or (BAP-tagged) SpoIIQ (XWB46) showed strong and specific reactivity to the antisera. Minor background

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
br>27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558