

SLC26A4 antibody

Cat#: orb417939 (CoA)

Product Name	SLC26A4 Antibody
Catalog Number	orb417939
Lot Number	G0914A
Immunogen	Recombinant Human Pendrin protein (592-655AA)
Host Species	Rabbit
Form	Liquid
Concentration	4.31mg/ml
Species Cross-reactivity	Human
Purification method	>95%,Protein G purified
Applications	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:200-1:500,IF:1:50-1:200
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage and Stability	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Conclusion	pass