

LEXSY *in vitro* Translation Cell Extract

Cat#: orb532694 (User Manual)

LEXSY *in vitro* Translation Cell Extract

For general laboratory use. Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Shelf Life: 12 months

Description:

The LEXSY cell extract for *in vitro* translation contains functional ribosomes and all essential components of the eukaryotic translation and folding machinery. For target mRNA generation the heterologous T7 RNA polymerase was added to the extracts.

Please note:

orb532694 is part of the LEXSY *in vitro* Translation Kit available as plasmid based or PCR based cell-free protein synthesis kit.