

**CRX antibody****Cat#: orb7145 (CoA)****Host:** Rabbit**Lot:** Q7690**Target Protein:** CRX1**IR:** Immunogen Range:51-150/299**Clonality:** Polyclonal**Isotype:** IgG**Entrez Gene:** [1406](#)**Swiss Prot:** [O43186](#)**Source:** KLH conjugated synthetic peptide derived from human CRX:51-150/299**Purification:** affinity purified by Protein A**Storage:** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.**Concentration:**1mg/ml**Applications:**

ELISA(1:5000-10000)

IHC-P(1:100-500)

IHC-F(1:100-500)

IF(1:100-500)

**Cross Reactive Species:**

Human

Mouse

Rat

Dog

Pig

Cow

Horse

Rabbit