



AID rabbit pAb

Cat#: orb767171 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	AID rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human AICDA. AA range:81-130
Specificity	AID Polyclonal Antibody detects endogenous levels of AID protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Activation-induced cytidine deaminase
Gene Name	AICDA
Cellular localization	Nucleus . Cytoplasm . Predominantly cytoplasmic (PubMed:21385873). In the presence of MCM3AP/GANP, relocalizes to the nucleus (By similarity)
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal



www.biorbyt.com

Concentration	1 mg/ml
Observed band	24kD
Human Gene ID	57379
Human Swiss-Prot Number	Q9GZX7
Alternative Names	AICDA; AID; Activation-induced cytidine deaminase; Cytidine aminohydrolase
Background	This gene encodes a RNA-editing deaminase that is a member of the cytidine deaminase family. The protein is involved in somatic hypermutation, gene conversion, and class-switch recombination of immunoglobulin genes. Defects in this gene are the cause of autosomal recessive hyper-IgM immunodeficiency syndrome type 2 (HIGM2). [provided by RefSeq, Feb 2009],



Western Blot analysis of Hela, HepG2 cells using AID Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-colon-cancer, antibody was diluted at 1:100