

ZBTB6 rabbit pAb**Cat#: orb773356 (Manual)**


For research use only. Not intended for diagnostic use.

Product Name	ZBTB6 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZBTB6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger and BTB domain-containing protein 6 (Zinc finger protein 482) (Zinc finger protein with interaction domain)
Gene Name	ZBTB6 ZID ZNF482
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	46kD
Human Gene ID	10773
Human Swiss-Prot Number	Q15916
Alternative Names	

Background

function: May be involved in transcriptional regulation., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 4 C2H2-type zinc fingers., tissue specificity: Widely expressed with highest levels in brain.,



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night