

Product Datasheet

Vazyme - FastPure FFPE DNA Isolation Kit (DC105-01)

Catalog Number	DC105-01
Category	Tools
Description	<p>FastPure FFPE DNA Isolation Kit is intended for isolating genomic DNA from formalin-fixed, paraffin-embedded (FFPE) tissue, overcoming the inhibitory effect caused by formalin crosslinking. The deparaffinization solution in this kit is safer than xylene. The DNA-specific spin column and the unique buffer system also ensure the high quality of obtained DNA. The whole extraction process takes only 20 min (digestion time not counted). The extracted genomic DNA is of good integrity, high purity, and stable quality and can be used for a wide range of downstream applications, including PCR, qPCR, and library preparation.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	<p>Proteinase K powder: Store at 2 ~ 8°C and adjust the shipping method according to the destination. Store dissolved Proteinase K at -30 ~ -15°C. Store the other components at 15 ~ 25°C and transport at room temperature.</p>
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)