

Product Datasheet

Vazyme - Magnetic Blood DNA Extraction Kit (Prepackaged) (DM102-01)

Catalog Number	DM102-01
Category	Tools
Description	<p>Magnetic Blood DNA Extraction Kit (Prepackaged) is suitable for the isolation and purification of high-quality genomic DNA from blood samples. The kit uses unique embedded superparamagnetic silicon-based magnetic beads. In a high-salt buffer system, nucleic acids instead of proteins and other impurities are adsorbed by hydrogen bonds and electrostatic binding. The magnetic beads that have adsorbed nucleic acids are washed to remove the remaining proteins and salts. When using low-salt buffer, nucleic acids are released from magnetic beads, so as to achieve the purpose of rapid separation and purification of nucleic acids. The kit is compatible with an automatic nucleic acid extraction instrument (Vazyme #VNP-32P) that leverages magnetic bead-based adsorption. Specially designed magnetic rods are used to adsorb, transfer, and release magnetic beads, allowing for automatic nucleic acid extraction and purification via the transfer of magnetic beads/nucleic acids. The DNA isolated with this kit are suitable for various downstream applications, including PCR, real-time PCR, library preparation, and biochip analysis, etc.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at 15 ~ 25°C and transport at room temperature. If ambient temperatures often exceed 25°C, we suggest storing Proteinase K at 2 ~ 8°C.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)