

Product Datasheet

Vazyme - VAHTS RNA Multiplex Oligos Set 2 for Illumina (N324-01)

Catalog Number	N324-01
Category	Tools
Description	VAHTS RNA Multiplex Oligos Set 1 for Illumina is specially designed for RNA library preparation for Illumina platforms. This kit contains VAHTA RNA Adapter-S for Illumina, 8 kinds of VAHTS i5 PCR Primers, and 12 kinds of VAHTS i7 PCR Primers. This kit provides 96 kinds of different dual-index combinations. Using this kit and VAHTS RNA Multiplex Oligos Set 2 for Illumina (Vazyme, #N324) together, can provide libraries with 384 kinds of different dual-index combinations.
Application notes	Quality Control: >99% purity confirmed by HPLC; nuclease-free components Reliable performance: Classical products for efficient ligation reactions
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15 °C and transport at ≤0 °C .
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Biorbyt Ltd.**

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)