

Product Datasheet

Vazyme - VAHTS mRNA Capture Beads 2.0 (N403-01)

Catalog Number	N403-01
Category	Tools
Description	VAHTS mRNA Capture Beads 2.0 are 1 µm Oligo(dT)-coupled paramagnetic microspheres intended for the isolation of poly(A) ⁺ RNA from purified total RNA. The magnetic separation technology allows for the isolation of intact mRNA in small sample volumes, eliminating the need for mRNA precipitation. The whole workflow can be completed in 1 h.
Application notes	High purity: Specific enrichment of mRNA, higher purity and reduced rRNA residue. Wide compatibility: Compatible with multiple species - eukaryotes such as animals, plants, fungi. Compatible with a wider range of input samples—10 ng-12.5 µg.
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at 2 ~ 8°C. Adjust the shipping method according to the destination.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)