

Product Datasheet

Vazyme - VAHTS Target Capture Universal Blockers and Post-PCR Primer Mix for Illumina-TS (NC101-01)

Description	VAHTS Target Capture Universal Blockers and Post-PCR Primer Mix for Illumina-TS contains Universal Blockers for Illumina-TS and Post-PCR Primer Mix for Illumina-TS. Universal Blockers for Illumina-TS are universal adapter blocking sequences specifically developed for common Illumina DNA libraries. They are applicable for single-indexed or dual-indexed DNA libraries (6 - 8 bp) and can prevent concatemerization of library molecules caused by annealing of adapters during hybridization, thus reducing non-specific hybridization and improving the on-target rate of hybridization capture. Post-PCR Primer Mix for Illumina-TS contains universal amplification primers specially developed for Illumina TruSeq DNA libraries. All the reagents provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of the experiment.
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Universal adapter blocking sequence developed for common DNA libraries on the Illumina platform. Suitable for single-end or paired-end 6 - 8 bp Barcode-labeled DNA libraries.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)