

Product Datasheet

Vazyme - Taq HS DNA Polymerase (P132-d2)

Catalog Number	P132-d2
Category	Tools
Description	Taq HS DNA polymerase is a hot-start Taq polymerase obtained by mixing Champagne Taq antibody with Taq DNA polymerase in an optimal ratio. Due to the unique thermo stability of Champagne Taq antibody, the activity of Taq HS DNA polymerase is still blocked at temperature up to 55°C, which minimizes non-specific amplification during the mixing and system heating. When the reaction is kept at 95°C for more than 30 sec, Champagne Taq antibody is completely inactivated and Taq enzyme activity is completely released, ensuring that the PCR system has extremely high amplification sensitivity and specificity. The activation of Taq HS DNA polymerase is not affected by pH, ionic strength, etc. It is applicable for various hot-start PCR and qPCR based on Taq DNA polymerase and can be used to amplify gene with low copy numbers from complex templates (genome and cDNA). It is the hot-start Taq enzyme of choice for PCR/qPCR molecular diagnostic reagents. Taq HS DNA polymerase has higher stability and detection rate.
Application notes	Hot-start technology: Superior specificity and sensitivity; convenient room-temperature reaction setup Good Stability: Stable performance after 30 freeze-thaw cycles or store at 4°C for 4 weeks Robust performance: Higher yield of target amplicons from low template amounts/complex templates
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)