

Product Datasheet

Vazyme - AccurSTART U+ One Step RT-qPCR Super PreMix (ONE TUBE) (Q621-03)

Catalog Number	Q621-03
Category	Tools
Description	<p>AccurSTART One Step RT-qPCR Super PreMix (ONE TUBE) is a one-tube probe-based RT-qPCR master mix suitable for single or multiple quantitative PCR detection using RNA as a template (such as RNA virus). The enzyme in this master mix uses a hot-start technology based on enzyme modification, with a new generation of high-affinity Taq enzyme antibodies, to achieve efficient blocking and precise release. Combined with one-step dedicated Reverse Transcriptase and optimized buffer, it supports pre-mixing of primers and probes in advance and the mixture can be stored for a long time at low temperature. The sample to be tested can be added directly when using, no additional tube opening/pipetting operation. This product is also compatible with fast programs, which greatly shortens the detection time.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and transport at ≤ 0°C.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)