

Product Datasheet

Vazyme - HiScript III RT SuperMix for qPCR (+gDNA wiper) (R323-00)

Catalog Number	R323-00
Category	Tools
Description	<p>HiScript III RT SuperMix for qPCR (+gDNA wiper) is an upgraded version of HiScript II Q RT SuperMix for qPCR (+gDNA wiper), including HiScript III Reverse Transcriptase, a new generation of reverse transcriptase with optimized Buffer. This kit further improves the efficiency of cDNA synthesis, and is suitable for two-step qRT-PCR detection. The 4 × gDNA wiper Mix in the kit can remove residual genomic DNA from the RNA template, ensuring more reliable quantitative results. It simplifies qPCR primer design without the need of spanning an exon-exon junction; 5 × HiScript III qRT SuperMix contains all the components required for the reverse transcription reaction. The reaction can be performed quickly by adding template RNA and RNase-free ddH₂O, and at the same time, the gDNA wiper function is terminated to ensure the integrity of the cDNA. It is compatible with dye-based and probe-based qPCR, enabling high-performance gene expression analysis.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15 °C and transport at ≤0 °C .
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)