

Product Datasheet

Vazyme - HiScript IV RT SuperMix for qPCR (+gDNA wiper) (R423-00)

Catalog Number	R423-00
Category	Tools
Description	HiScript IV RT SuperMix for qPCR (+gDNA wiper) is an upgraded version of HiScript III RT SuperMix for qPCR (+gDNA wiper), including a new generation of HiScript IV Reverse Transcriptase and Buffer optimized for reverse transcription. This kit further improves the synthesis efficiency of cDNA, making it a better choice for reverse transcription of low-input, low-expression or degraded RNA templates. 5 × gDNA wiper Mix can completely remove the genomic contamination in the RNA template, so there is no need to design intron-spanning qPCR primers. 4 × HiScript IV qRT SuperMix contains all the components required for the reverse transcription, just add template RNA and RNase-free ddH₂O to start the reaction.
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15 °C and transport at ≤0 °C .
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)