

Product Datasheet

Vazyme - FreeZol Reagent (R711-00)

Catalog Number	R711-00
Category	Tools
Description	<p>FreeZol Reagent is widely applicable to the extraction of total RNA from cultured cells, animal tissues, and simple plant tissues. Compared with the conventional Trizol method, this product features a simple procedure that can be performed at room temperature with no need to use chloroform for phase separation. In addition, this product ensures the integrity and purity of the extracted RNA by precipitating proteins, DNA, polysaccharides and other impurities in the organic phase while retaining RNA in the upper aqueous phase. The whole procedure can be completed in 1 h. The obtained total RNA can be used directly for RT-PCR, qRT-PCR, Northern blot, Dot blot, in vitro translation, NGS, and other molecular biology experiments.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at 2~8°C and transport them at room temperature.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)