

Product Datasheet

Vazyme - VeZol-Pure Total RNA Isolation Kit (RC202-01)

Catalog Number	RC202-01
Category	Tools
Description	<p>VeZol-Pure Total RNA Isolation Kit is applicable for isolating high-quality Total RNA from various samples such as cultured cells, animal and plant tissues. This kit combines phenol/guanidine-based lysis of samples with silica membrane-based purification of total RNA, with the advantages of strong lysing ability, wide compatibility, fast, high yields, and good purity. Samples are fully lysed with VeZol Reagent, and RNA is specifically adsorbed from the aqueous phase by silica membrane. The obtained total RNA can be directly used for molecular biology experiments such as RT-PCR, qRT-PCR, Northern Blot, Dot Blot, in vitro translation, and high-throughput sequencing.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store BOX 1 at 2 ~ 8°C and protect from light. Transport at room temperature. Store BOX 2 at 15 ~ 25°C and transport at room temperature.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)