

Product Datasheet

Vazyme - FastPure Viral DNA/RNA Mini Kit Pro (RC323-00)

Catalog Number	RC323-00
Category	Tools
Description	<p>FastPure Viral DNA/RNA Mini Kit Pro is a column extraction kit, which can extract viral DNA/RNA from multiple types of samples. The kit uses a unique lysis system, without high temperature incubation, and can significantly improve the recovery rate of trace nucleic acid with carrier RNA. After specific adsorption on silica gel membrane, viral DNA/RNA can be quickly and efficiently purified. For a single sample, it only takes 12 minutes to finish extraction. And the products can be directly used for downstream related experiments such as reverse transcription, PCR, Real-Time PCR, NGS and Northern blot.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Box 1: Store at -30 ~ -15°C and transport at $\leq 0^{\circ}\text{C}$. Box 2: Store at 15 ~ 25°C and transport at room temperature.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)