

## Product Datasheet

### Vazyme - LAMP Fluorescent Dye (RP001-01)

#### Description

LAMP Fluorescent Dye is an intercalating dye with the strong fluorescent signal and little inhibition of the reaction. It is applicable for real-time fluorescent detection of LAMP reactions. This product is a 50 × solution with DMSO as the solvent, compatible with a variety of fluorescent quantitative PCR instruments. It is recommended that the final concentration of the dye be 0.1 - 1×, and please select the SYBR/FAM detection channel when using it.

#### Storage

Store at -30 ~ -15°C and protect from light. Transport at ≤0°C. ▲ Avoid repeated freezing and thawing.

#### Note

For research use only

#### Application notes

Robust performance: Strong fluorescent signal, little inhibition on reaction  
Optimized formulation: Excellent amplification curve  
Broad Compatible: Applicable to most instruments

#### Expiration Date

12 months from date of receipt.



#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)