

Product Datasheet

Type III Collagen Antibody (orb1251420)


Description	Type III Collagen Antibody
Species/Host	Goat
Reactivity	Bovine, Canine, Equine, Gallus, Guinea pig, Hamster, Human, Mammal, Mouse, Porcine, Primate, Rabbit, Rat, Rodent
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Type III collagen
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.4 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Type III Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, flow cytometry, and western blotting assays.
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Human gastric cancer tissue was stained with Goat Anti-Type III Collagen-UNLB (Cat. No. orb1251420) followed by Swine Anti-Goat IgG (H + L), Human/Rat/Mouse SP ads-HRP.



Frozen osteopontin-null mouse muscle section was stained with Goat Anti-Type III Collagen-UNLB (Cat. No. orb1251420) followed by a secondary antibody.




Rat kidney section was stained with Goat Anti-Type III Collagen-UNLB (Cat. No. orb1251420) followed by a secondary antibody and DAB and hematoxylin.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Frozen kidney section following unilateral ureteral obstruction was stained with Goat Anti-Type III Collagen-UNLB (Cat. No. orb1251420) followed by a secondary antibody.



Paraffin embedded doxorubicin treated rat kidney section was stained with Goat Anti-Type III Collagen-UNLB (Cat. No. orb1251420) followed by a secondary antibody and DAB.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)