

Product Datasheet

SYBR Color qPCR SuperMix (Low ROX Premixed) (orb1736810)

Catalog Number orb1736810

Category Tools

Description 1) Quickly get results, saving up to 50% of the time ||2) Optimized ready-to-use pre mix for rapid PCR reactions ||3) Accurate detection of various starting amounts of templates, stable amplification, quantitative results with high repeatability ||4) Balanced K⁺ and NH₄⁺ ion ratios, ensure high sensitivity and high specificity ||5) Track the pipetting process and reduce pipetting errors through the color change reaction between dyes. . SYBR Color qPCR SuperMix is an optimized 2x real-time PCR master mix containing HotStarTaq DNA Polymerase, SYBR Green fluorescent dye, dNTP and Mg²⁺, low Rox. In addition, the balanced K⁺ and NH₄⁺ ion ratios in the buffer promote specific primer annealing. To ensure a highly sensitive and specific PCR reaction, the reaction can be initiated by simply adding the primer and cDNA template to the ready-to-use PCR master mix. Greatly simplify the operation process and reduce the probability of pollution. The unique PCR buffer ensures sensitive qPCR on all Real-Time PCR instruments without optimization. |In addition, the product realizes the visualization of the sample addition process through the discoloration reaction after the template is added, which greatly improves the sample addition efficiency and avoids sample addition errors

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)