

Product Datasheet

50×ROX Reference Dye (orb1736814)

Catalog Number	orb1736814
Category	Tools
Description	ROX Reference Dye is generally used on ABI, Stratagene and other companies' Real Time PCR amplification machines to adjust the difference between PCR tubes and tubes caused by PCR sample loading errors. This product is suitable for Applied Biosystems 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast; Applied Biosystems StepOne, StepOnePlus and other fluorescent quantitative PCR instruments using high concentration ROX. ROX Reference Dye is dissolved in a buffer system of 20 mM Tris-HCl (pH 8.4), 0.1 mM EDTA, and 0.01% Tween-20.. .
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)