

Product Datasheet

KBA.62 (Melanoma Associated Antigen) (orb1945630)

Description	KBA.62, also known as Melanoma Associated Antigen, is used to detect an antigen present in melanocytic tumours, such as melanomas, due to its proven sensitivity and specificity. The antibody can also be used to distinguish between junctional nevus and intradermal nevus cells, and fetal melanocytes versus normal adult melanocytes. Studies have shown KBA.62 to be highly useful in differentiating between metastatic amelanotic melanoma and a number of poorly differentiated carcinomas, large cell lymphomas, sarcomas, and spindle cell carcinomas.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Storage	Store at 2-8°C. To ensure stability, immediately replace vial back in the refrigerator after each use
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	IHC062
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)