

Product Datasheet

Rat 25 Hydroxy Vitamin D (25-OH-D) ELISA Kit (orb2307978)

Catalog Number orb2307978

Description This assay employs the competitive inhibition enzyme immunoassay technique.

The microtiter plate provided in this kit has been pre-coated with Rat 25-OH-D protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rat 25-OH-D. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm.

The concentration of Rat 25-OH-D in the samples is then determined by

comparing the OD of the samples to the standard curve.

Reactivity Rat

Range 3.13-200 ng/mL

Concentration 200 ng/mL

Note For research use only

Application notes standard: 200 ng/mL. Test principle: This assay employs the competitive

inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Rat 25-OH-D protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rat 25-OH-D. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the

addition of sulphuric acid solution and the color change is measured

spectrophotometrically at a wavelength of $450 \text{nm} \pm 10 \text{nm}$. The concentration of Rat 25-OH-D in the samples is then determined by comparing the OD of the

samples to the standard curve

Sample Types serum, plasma and other biological fluids

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558





Assay Time 2.5h

Sensitivity 0.93 ng/mL

Expiration Date Please enquire.

 $\begin{aligned} & \text{Email: } \underline{\textbf{info@biorbyt.com}}, \ \underline{\textbf{support@biorbyt.com}} \\ & \text{Phone: } \underline{\textbf{+1 (415) 906-5211}} \mid \text{Fax: } \underline{\textbf{+1 (415) 651-8558}} \end{aligned}$