

Product Datasheet

Horse 5-Hydroxytryptamine (5-HT) ELISA Kit (orb2308024)

Catalog Number orb2308024

Description This assay employs the competitive inhibition enzyme immunoassay technique.

The microtiter plate provided in this kit has been pre-coated with Horse 5-HT protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Horse 5-HT. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of $450 \, \mathrm{nm} \pm 10 \, \mathrm{nm}$.

The concentration of Horse 5-HT in the samples is then determined by

comparing the OD of the samples to the standard curve.

Reactivity Equine

Range 14.07-900 ng/mL

Concentration 900 ng/mL

Note For research use only

Application notes standard: 900 ng/mL. Test principle: This assay employs the competitive

inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Horse 5-HT protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Horse 5-HT. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the

addition of sulphuric acid solution and the color change is measured

spectrophotometrically at a wavelength of $450 \text{nm} \pm 10 \text{nm}$. The concentration of Horse 5-HT in the samples is then determined by comparing the OD of the

samples to the standard curve

Sample Types serum, plasma, tissue homogenates





Assay Time 2.5h

Sensitivity 5.85 ng/mL

Expiration Date Please enquire.

68 TW Alexander Drive, Durham, NC, 27713, United States

$$\label{eq:phone:mail:model} \begin{split} &\text{Email: } \underline{info@biorbyt.com}, \ \underline{support@biorbyt.com} \\ &\text{Phone: } \underline{+1 \ (415) \ 906-5211} \ \big| \ \text{Fax: } \underline{+1 \ (415) \ 651-8558} \end{split}$$