

## Product Datasheet

### Direct Extraction Buffer (orb2646461)

<b>Catalog Number</b>	orb2646461
<b>Category</b>	Tools
<b>Description</b>	Direct Extraction Buffer allows an easy and fast extraction of DNA and RNA directly from blood, swabs and animal- or plant tissue. The buffer is optimized for use in combination with Direct PCR or RT-PCR master mixes like qPCR ProbesMaster (PCR-396) or SCRIPT Direct RT-qPCR ProbesMaster (PCR-528).\n\nThe mix allows DNA and RNA preparation within 3-5 minutes and with a minimum of pipetting steps. It is especially recommended for:
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	5 x
<b>Storage</b>	store at 4°C or -20°C
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)